

Sound Pro™ Combination Table – Ultrasound Table

058-710

058-715

BIODEX SOUND PRO COMBINATION TABLE, 115 VAC

BIODEX SOUND PRO COMBINATION TABLE, 230 VAC

- Dimensions:
 - Overall: 71" l x 30" w (180.3 x 76.2 cm); 35" w (88.9 cm) with side rails
 - Cardiac Scanning Cutout: 9.5" w x 15.75" l (24.1 x 40 cm)
Drop-down cushion releases from either side of the table.
 - Sonographer Cutout: 8" w x 14.75" l (20.3 x 37.4 cm)
Adjustable sonographer cushion drops down or folds up for maximum access to the patient.
 - Drop-Down Leg Section: 12.6" l (32 cm)
- Motions:
 - Height Adjustable: 23" to 39" (58.4 to 99 cm)
 - Trendelenburg: 0° to ±15°
 - Fowler Back: 0° to 80° infinitely adjustable
- Controls:
 - Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
 - Foot Control: Activates height, Trendelenburg motions and Fowler positioning
- Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system
- Tabletop: Three primary sections: torso section features sonographer and patient cutouts with Fowler positioning; center section remains fixed and leg section drops down to 40° and 80° for stirrup access
 - Drop-Down Cardiac Scanning Cushion: Release from either side of the table
 - Adjustable Sonographer Cushion: drops down or folds up for maximum access to the patient
- Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Upholstery Color: Graphite
- Patient Restraints: One body strap
- Finish: Powder coat
- Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
- Weight: 350 lb (159 kg)
- Power: 115 VAC or 230 VAC
- Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN 60601-1 standards.
- Warranty: Two years parts and labor



Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, NY, USA, 11967-4704, Phone: 001-631-924-9000, Email: info@biodex.com, www.biodex.com