

## Vasc Pro™ Vascular Ultrasound Table – For Long-Leg Vascular Studies

058-732

058-733

**BIODEX VASC PRO VASCULAR ULTRASOUND TABLE, 115 VAC**

**BIODEX VASC PRO VASCULAR ULTRASOUND TABLE, 230 VAC**

- Dimensions:
  - Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails
- Motions:
  - Height Adjustable: 23" to 39" (58.4 to 99 cm)
  - Trendelenburg: 0° to 15°
  - Reverse Trendelenburg: 0° to 30°
  - Fowler Back: 0° to 80° infinitely adjustable
- Controls:
  - Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
  - Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning
- Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system
- Tabletop: Two primary sections: torso section with Fowler positioning; lower body section remains fixed
- Foot Support: Adjustable every 2", up to 14"
- Mattress: Torso section 2" (5 cm) thick, lower body section 3" (7.6 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Upholstery Color: Graphite
- Patient Restraints: Three body straps
- Finish: Powder coat
- Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
- Weight: 335 lb (152 kg)
- Power: 115 VAC or 230 VAC
- Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1.M90 and EN60601-1 standards.
- Warranty: Two years parts and labor



**Biodex Medical Systems, Inc.**

20 Ramsey Road, Shirley, NY, USA, 11967-4704, Phone: 001-631-924-9000, Email: [info@biodex.com](mailto:info@biodex.com), [www.biodex.com](http://www.biodex.com)