

Brachytherapy C-Arm Table - 810 – For Seed Implantation Procedures

058-810

BIODEX BRACHYTHERAPY C-ARM TABLE – 810, 115 VAC

058-815

BIODEX BRACHYTHERAPY C-ARM TABLE – 810, 230 VAC

- Dimensions: 84.5" l x 25.5" w (214.6 x 64.8 cm) with OR accessory rails
- Tabletop: 53" l x 21" w (134.6 x 53.3 cm); 84" l (213.4 cm) with extension
- Lightweight Extension (for patient transfer only, not for imaging): 26.5" l x 21" w (67.3 x 53.3 cm)
- Radiolucent Extension: 31" l x 25.5" w (78.7 x 64.8 cm) with OR accessory rails
- Radiolucent Area: 54" l x 17.5" w (137.2 x 44.5 cm)
- Tabletop Material: Three Fibresin™ sections – each with individual 2" mattress
- Mattress: Individual, seamless, 2" thick (5 cm), one for each section
- Motions:
 - Height Adjustable: 34" to 43.5" (86.4 to 110.5 cm)
 - Tabletop X Motion (head-to-toe): 20" (51 cm)
 - Tabletop Y Motion (side-to-side): 10" (25.4 cm)
 - Trendelenburg: 0° to 20°
 - Reverse Trendelenburg: 0° to 20°
 - Lateral Roll: 0° to 20°
 - Fowler Positioning: Adjusts from 0° to 80° with touch handle lock/unlock
- Controls:
 - Hand Control: activates all motions
 - Foot Control: activates all motions
- Optional:
 - X-Ray Film Tray: Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table.
 - X-Ray Grid: Removable X-Ray grid, 103 line, 8:1 ratio
- Accessory Rails: Standard OR accessory rails 1.125" x .375" (2.86 x .95 cm) on table (rails also on table imaging extension)
- Patient Restraints: Two body straps, two ankle straps
- Finish: Powder coat
- Patient Capacity:
 - In urologic position: 500 lb (227 kg)
 - In use with radiolucent extension: 400 lb (181.4 kg)
 - Weight tested to four times the patient load rating
- Shipping Weight: 877 lb (398.6 kg)
- Power: 115 VAC or 230 VAC and battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.
- Certifications: ETL and cETL listed to UL 60601-1®, CAN/CSA© C22.2 No.: 601.1.M90, EN 60601-1, IEC 60601-2-46 and CE marked.
- CDRH registered and FDA Device Listed
- Warranty: Two years parts and labor



Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, NY, USA, 11967-4704, Phone: 001-631-924-9000, Email: info@biodex.com, www.biodex.com