

NEW Pain Management C-Arm Table - 870 – Ideal for Pain Care Applications

058-870

BIODEX PAIN MANAGEMENT C-ARM TABLE – 870 Contoured Top, 115 VAC

058-870-10

BIODEX PAIN MANAGEMENT C-ARM TABLE – 870-10 Rectangular Top, 115 VAC

058-875

BIODEX PAIN MANAGEMENT C-ARM TABLE – 870 Contoured Top, 230 VAC

058-875-10

BIODEX PAIN MANAGEMENT C-ARM TABLE – 870-10 Rectangular Top, 230 VAC

- Dimensions:
Overall: 97" l x 26" w (246 x 66 cm) with OR accessory rails
Tabletop: 97" l x 24" w (246 x 61 cm) – choice of classic contoured or rectangular design
- Radiolucent Area: 70" l x 24" w (178 x 61 cm)
- Tabletop Material: Carbon fiber with integral head section
- Mattress: Seamless, 2" thick (5 cm)
- Attenuation: 1.2 mm Aluminum equivalence
- Motions:
Height Adjustable: 29.5" to 39.5" (75 to 100 cm)
Trendelenburg: 0° to 20°
Reverse Trendelenburg: 0° to 20°
Isocentric Lateral Roll: 15° clockwise or counterclockwise
- Controls:
Hand Control: activates height, isocentric lateral roll and Trendelenburg motions
Foot Control (optional): activates height, lateral roll and Trendelenburg motions
- Accessory Rails: Standard OR accessory rails 30" x 1.12" x .375" (75 x 2.86 x .95 cm) mounted near foot end of table
- Patient Restraints: Three body straps
- Finish: Stainless Steel
- Patient Capacity: 500 lb (227 kg)
- Shipping Weight: 749 lb (340 kg)
- Power: 115 VAC or 230 VAC
- Certifications: ANSI/AAMI ES60601-1:2005 +C1 (2009) +A1 (2012) +A2 (2010); CAN/CSA C22.2 No: 60601-1:14; CAN/CSA-C22.2 NO. 60601-2-46:12; IEC 60601-1-2-46:2016.
- CDRH registered and FDA Device Listed
- Warranty: Two years parts and labor



Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, NY, USA, 11967-4704, Phone: 001-631-924-9000, Email: info@biodex.com, www.biodex.com